



MELTING POINT TESTER

SSP

## The ERWEKA SSP measures the melting point of suppository samples.

It consists of a graduated tube with an integrated test chamber made of glass. The sample to be tested is placed in the spiral-shaped glass test basket inside the test chamber, which is surrounded by a water jacket. The water inside the jacket can easily be heated with the optional EST 2 circulation heater.



## OPTIONS

- | ERWEKA circulation heater EST 2
- | IQ / OQ documents

### MELTING POINT TESTER SSP

## TECHNICAL DATA

---

**Weight**

5 kg

---

**Dimensions (H x W x D)**

750 x 240 x 280 mm

---

**Test stations**

1

---

**Temperature measuring range**

25 - 50 °C ± 0.2 °C

---

**Heating**

Thermostat of EST 2 / ED-5 / Corio CD-B5 type

---